




Recombinant Mouse Interleukin-36 beta protein (Il36b) (Active)

Product Code	CSB-AP003111MO
Abbreviation	Recombinant Mouse Il36b protein (Active)
Uniprot No.	Q9D6Z6
Form	Lyophilized powder
Storage Buffer	Lyophilized from a 0.2 µm filtered 20 mM Tris, 300 mM NaCl, pH 8.0, 5 % trehalose
Product Type	Interleukin
Immunogen Species	Mus musculus (Mouse)
Biological Activity	Fully biologically active when compared to standard. The specific activity determined by its ability in a functional ELISA. Immobilized rMuIL-36β at 1 ug/mL can bind recombinant murine IL-1 Rrp2 with a range of 0.15-5 ug/mL.
Purity	≥95% as determined by SDS-PAGE.
Sequence	MMAFPPQSCV HVLPPKSIQM WEPNHNTMHG SSQSPRNYRV HDSQQMVVWL TGNTLTAVPA SNNVKPVILS LIACRDTEFQ DVKKGNLVFL GIKNRNLCFC CVEMEGKPTL QLKEVDIMNL YKERKAQKAF LFYHGIEGST SVFQSVLYPG WFIATSSIER QTIILTHQRG KLVNTNFYIE SEK
Research Area	Immunology
Source	E.coli
Target Names	Il36b
Expression Region	1-183aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag-Free
Mol. Weight	20.9 kDa
Protein Length	Full Length
PubMed ID	11991723; 16141072; 19717513; 21860022; 21965679
Image	


Endotoxin

Less than 1.0 EU/μg as determined by LAL method.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.