




# Recombinant Human Stromal cell-derived factor 1 protein (CXCL12) (Active)

<b>Product Code</b>	CSB-AP000741HU
<b>Abbreviation</b>	Recombinant Human CXCL12 protein, partial (Active)
<b>Uniprot No.</b>	P48061
<b>Storage Buffer</b>	20 mM PB pH 7.0, 130 mM NaCl
<b>Product Type</b>	Chemokines
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Biological Activity</b>	Fully biologically active when compared to standard. The biological activity determined by a chemotaxis bioassay using PHA and rHuIL-2 activated human peripheral blood T-lymphocytes is in a concentration range of 20-80 ng/ml.
<b>Purity</b>	≥97% as determined by SDS-PAGE.
<b>Sequence</b>	KPVLSYRCP CRFFESHVAR ANVKHLKILN TPNCALQIVA RLKNNNRQVC IDPKLKWIQE YLEKALNK
<b>Research Area</b>	Immunology
<b>Source</b>	E.Coli
<b>Target Names</b>	CXCL12
<b>Expression Region</b>	22-89aa
<b>Tag Info</b>	Tag-Free
<b>Mol. Weight</b>	8.0 kDa
<b>Protein Length</b>	Full Length of Mature Protein of Isoform Alpha
<b>PubMed ID</b>	7490086; 16626895; 14702039; 16344560; 15164054; 15489334; 8752281; 9427609; 10446158; 11069075; 11859124; 14525775; 16107333; 15741341; 16725153; 18802065; 19255243; 22457824; 9384579; 9618518; 10954912; 15588815; 17264079; 17357154; 18799424; 19551879;

## Image



**Endotoxin**

Less than 1.0 EU/μg as determined by LAL method.

**Shelf Life**

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.