




Recombinant Macaca mulatta C-X-C motif chemokine 10 protein (CXCL10)

Product Code	CSB-AP001031MOW
Abbreviation	Recombinant Rhesus macaque CXCL10 protein (Active)
Uniprot No.	Q8MIZ1
Storage Buffer	0.2 μm filtered PBS, pH 7.4 ,lyophilized
Product Type	Chemokines
Immunogen Species	Macaca mulatta (Rhesus macaque)
Biological Activity	Fully biologically active when compared to standard. The biological activity determined by a chemotaxis bioassay using human peripheral blood T-lymphocytes is in a concentration range of 10- 50 ng/ml.
Purity	≥97% as determined by SDS-PAGE.
Sequence	IPLSRTVRCT CISISNQPVN PRSLEKLEII PPSQFCPHVE IIATMKKKGE KRCLNPESKA IKNLLKAVSK ERSKRSP
Research Area	Immunology
Source	E.Coli
Target Names	CXCL10
Expression Region	22-98aa
Tag Info	Tag-Free
Mol. Weight	8.7 kDa
Protein Length	Full Length of Mature Protein

Image



Endotoxin	Less than 1.0 EU/μg as determined by LAL method.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.