


Recombinant Rhesus Macaque Eotaxin protein (CCL11)

Product Code	CSB-AP001041MOW
Abbreviation	Recombinant Rhesus macaque CCL11 protein (Active)
Uniprot No.	Q8MIT7
Storage Buffer	Lyophilized from a 0.2 μ m filtered concentrated solution in 2 $\times$ PBS, pH 6.9.
Product Type	Chemokines
Immunogen Species	Macaca mulatta (Rhesus macaque)
Biological Activity	Fully biologically active when compared to standard. The biological activity determined by a chemotaxis bioassay using human peripheral blood eosinophils is in a concentration range of 0.1-10.0 ng/ml.
Purity	$\geq$ 98% as determined by SDS-PAGE.
Sequence	GPDSVATTCC FTLTNKKIPL QRLESYRRII SGKCPQKAVI FKTKLAKDIC ADPKKKWVQD SMKYLDRKSP TPKP
Research Area	Immunology
Source	E.Coli
Target Names	CCL11
Expression Region	24-97aa
Tag Info	Tag-Free
Mol. Weight	8.4 kDa
Protein Length	Full Length of Mature Protein

Image



Endotoxin	Less than 1.0 EU/ μ g as determined by LAL method.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.