




Recombinant Mouse C-C motif chemokine 19 protein (Ccl19) (Active)

Product Code	CSB-AP001321MO
Abbreviation	Recombinant Mouse Ccl19 protein (Active)
Uniprot No.	O70460
Form	Lyophilized powder
Storage Buffer	Lyophilized from a 0.2 µm filtered PBS, pH 7.4
Product Type	Chemokine
Immunogen Species	Mus musculus (Mouse)
Biological Activity	Fully biologically active when compared to standard. The biological activity determined by a chemotaxis bioassay using human mature dendritic cells is in a concentration range of 10-100 ng/ml.
Purity	>97% as determined by SDS-PAGE.
Sequence	GANDAEDCCL SVTQRPIPGN IVKAFRYLLN EDGCRVPAVV FTTLRGYQLC APPDQPWVDR IIRRLKKSSA KNKGNSTRRS PVS
Research Area	Immunology
Source	E.coli
Target Names	Ccl19
Expression Region	26-108aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag-Free
Mol. Weight	9.2 kDa
Protein Length	Full Length of Mature Protein
PubMed ID	9653094; 15489334

Image



**Endotoxin**

Less than 1.0 EU/ μ g as determined by LAL method.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.