




Recombinant Mouse C-C motif chemokine 22 protein (Ccl22)

Product Code	CSB-AP001351MO
Abbreviation	Recombinant Mouse Ccl22 protein
Uniprot No.	O88430
Storage Buffer	0.2 μm filtered concentrated solution in 20 mM PB, pH 7.4, 150 mM NaCl, lyophilized
Product Type	Chemokines
Immunogen Species	Mus musculus (Mouse)
Biological Activity	Not Test
Purity	>97% as determined by SDS-PAGE.
Sequence	GPYGANVEDS ICCQDYIRHP LPSRLVKEFF WTSKSCRKPG VVLITVKNRD ICADPRQVWV KKLLHKLS
Research Area	Immunology
Source	E.Coli
Target Names	Ccl22
Expression Region	25-92aa
Tag Info	Tag-Free
Mol. Weight	7.8 kDa
Protein Length	Full Length of Mature Protein
PubMed ID	9687523; 15489334

Image



Endotoxin	Less than 1.0 EU/μg as determined by LAL method.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.