




Recombinant Rat BRAK protein (Cxcl14)

Product Code	CSB-AP001481RA
Abbreviation	Recombinant Rat Cxcl14 protein (Active)
Uniprot No.	Q8K453
Storage Buffer	Lyophilized from a 0.2 μm filtered concentrated solution in 30 % Acetonitrile and 0.1 % TFA.
Product Type	Chemokines
Immunogen Species	Rattus norvegicus (Rat)
Biological Activity	Fully biologically active when compared to standard. The biological activity determined by a chemotaxis bioassay using human monocytes is in a concentration range of 1.0-10 ng/ml.
Purity	≥95% as determined by SDS-PAGE.
Sequence	SKCKCSRKGP KIRYSDVKKL EMKPKYPHCE EKMVIITTKS MSRYRGQEHC LHPKLQSTKR FIKWYNWNE KRRVYEE
Research Area	Immunology
Source	E.Coli
Target Names	Cxcl14
Expression Region	23-99aa
Tag Info	Tag-Free
Mol. Weight	9.4 kDa
Protein Length	Full Length of Mature Protein
PubMed ID	15057822

Image



Endotoxin	Less than 1.0 EU/μg as determined by LAL method.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.