




Recombinant Human Insulin-like growth factor-binding protein 3 protein (IGFBP3) (Active)

Product Code	CSB-AP002491HU
Abbreviation	Recombinant Human IGFBP3 protein (Active)
Uniprot No.	P17936
Form	Lyophilized powder
Storage Buffer	Lyophilized from a 0.2 µm filtered PBS, pH 7.4
Product Type	Growth Factor
Immunogen Species	Homo sapiens (Human)
Biological Activity	Fully biologically active when compared to standard. The ED50 as determined by inhibiting IGF-II induced proliferation of serum free human MCF-7 cells is less than 200 ng/ml, corresponding to a specific activity of $>5.0 \times 10^3$ IU/mg in the presence of 15 ng/ml of rHuIGF-II.
Purity	≥98% as determined by SDS-PAGE.
Sequence	GASSAGLGPVVRCEPCDARALAQCAPPPAVCAELVREPGCGCCLTCALSEG QPCGIYTERCGSGLRCQPSPDEARPLQALLDGRGLCVNASAVSRLRAYLLPA PPAPGNASESEEDRSAGSVESPSVSSTHRVSDPKFHPLHSKIIKKGHAKDSQ RYKVDYESQSTDTQNFSSSESKRETEYGPCRREMEDTLNHLKFLNVLSPRGVHI PNCDDKKGFYKKKQCRPSKGRKRGFCWCVDKYGQPLPGYTTKGKEDVHCYS MQSK
Research Area	Cancer
Source	E.coli
Target Names	IGFBP3
Expression Region	28-291aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag-Free
Mol. Weight	28.8 kDa
Protein Length	Full Length of Mature Protein
PubMed ID	2464130; 1695633; 7508771; 12853948; 12690205; 15489334; 1697583; 1726837; 15340161; 14561895; 14760718; 16335952; 19159218; 19139490; 20353938; 26091039; 16959974

Image


Endotoxin

Less than 1.0 EU/μg as determined by LAL method.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.