




# Recombinant Mouse Tumor necrosis factor protein (Tnf), partial (Active)

<b>Product Code</b>	CSB-AP003471MO
<b>Abbreviation</b>	Recombinant Mouse Tnf protein, partial (Active)
<b>Uniprot No.</b>	P06804
<b>Form</b>	Lyophilized powder
<b>Storage Buffer</b>	Lyophilized from a 0.2 µm filtered PBS, pH 7.2
<b>Product Type</b>	Tumor Necrosis Factor
<b>Immunogen Species</b>	Mus musculus (Mouse)
<b>Biological Activity</b>	Fully biologically active when compared to standard. The ED50 as determined by a cytotoxicity assay using murine L929 cells is less than 0.1 ng/ml, corresponding to a specific activity of $>1.0 \times 10^7$ IU/mg in the presence of actinomycin D.
<b>Purity</b>	>98% as determined by SDS-PAGE.
<b>Sequence</b>	M+LRSSSQNSS DKPVAHVVAN HQVEEQLEWL SQRANALLAN GMDLKDNQLV VPADGLYLVY SQVLFKGGQC PDYVLLTHTV SRFAISYQEK VNLLSAVKSP CPKDTPEGAE LKPWYEPIYL GGVFQLEKGD QLSAEVNLPK YLDFAESGQV YFGVIAL
<b>Research Area</b>	Cancer
<b>Source</b>	E.coli
<b>Target Names</b>	Tnf
<b>Expression Region</b>	80-235aa
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag-Free
<b>Mol. Weight</b>	17.4 kDa
<b>Protein Length</b>	Partial
<b>PubMed ID</b>	2836146; 3898078; 2419912; 2989794; 3040015; 3684584; 7560085; 9089109; 14656967; 2777790; 2268312; 3349526; 10089307

Image


**Endotoxin**

Less than 1.0 EU/μg as determined by LAL method.

**Reconstitution**

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

**Shelf Life**

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.