




# Recombinant Mouse Uteroglobin protein (Scgb1a1) (Active)

<b>Product Code</b>	CSB-AP000561MO
<b>Abbreviation</b>	Recombinant Mouse Scgb1a1 protein (Active)
<b>Uniprot No.</b>	Q06318
<b>Form</b>	Lyophilized powder
<b>Storage Buffer</b>	Lyophilized from a 0.2 µm filtered 20 mM PB, pH 7.4, 150mM NaCl
<b>Product Type</b>	Other
<b>Immunogen Species</b>	Mus musculus (Mouse)
<b>Biological Activity</b>	Fully biologically active when compared to standard. The ED50 as determined by the ability of the immobilized protein to support the adhesion of the A549 human lung carcinoma cells is less than 5.0 µg/ml, corresponding to a specific activity of >200 IU/m.
<b>Purity</b>	≥98% as determined by SDS-PAGE.
<b>Sequence</b>	DICPGFLQVL EALLMESESG YVASLKPFNP GSDLQNAGTQ LKRLVDTLPQ ETRINIMKLT EKILTSPLCK QDLRF
<b>Research Area</b>	Immunology
<b>Source</b>	E.Coli
<b>Target Names</b>	Scgb1a1
<b>Expression Region</b>	22-96aa
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag-Free
<b>Mol. Weight</b>	8.4 kDa
<b>Protein Length</b>	Full Length of Mature Protein
<b>PubMed ID</b>	8398159; 8020953; 8440203; 7916613; 15489334; 21183079

## Image





---

**Endotoxin**

Less than 1.0 EU/ $\mu$ g as determined by LAL method.

---

**Reconstitution**

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

---

**Shelf Life**

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.