




Recombinant Human Fms-related tyrosine kinase 3 ligand protein (FLT3LG), partial (Active)

Product Code	CSB-AP002201HU
Abbreviation	Recombinant Human FLT3LG protein, partial (Active)
Uniprot No.	P49771
Storage Buffer	0.2 μm filtered PBS, pH 7.0 ,lyophilized
Product Type	Others
Immunogen Species	Homo sapiens (Human)
Biological Activity	Fully biologically active when compared to standard. The ED50 as determined by a cell proliferation assay using human AML5 cells is less than 1.0 ng/ml, corresponding to a specific activity of $>1.0 \times 10^6$ IU/mg.
Purity	$>97\%$ as determined by SDS-PAGE.
Sequence	TQDCSFQHSP ISSDFAVKIR ELSDYLLQDY PVTVASNLQD EELCGGLWRL VLAQRWMERL KTVAGSKMQG LLERVNTEIH FVTKCAFQPP PSCLRFVQTN ISRLQETSE QLVALKPWIT RQNFSCRLEL QCQPDSSTLP PPWSPRPLEA TAPTA
Research Area	Immunology
Source	E.Coli
Target Names	FLT3LG
Expression Region	27-181aa
Tag Info	Tag-Free
Mol. Weight	17.6 kDa
Protein Length	Partial
PubMed ID	8145851; 8180375; 7566977; 14702039; 15057824; 15489334; 10881197

Image



Endotoxin

Less than 1.0 EU/μg as determined by LAL method.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.