



S1PR1 Antibody, FITC conjugated

Product Code	CSB-PA020650LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P21453
Immunogen	Recombinant Human Sphingosine 1-phosphate receptor 1 protein (311-382AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Sphingosine 1-phosphate receptor 1 (S1P receptor 1) (S1P1) (Endothelial differentiation G-protein coupled receptor 1) (Sphingosine 1-phosphate receptor Edg-1) (S1P receptor Edg-1) (CD antigen CD363), S1PR1, CHEDG1 EDG1
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	S1PR1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.