



ATP2A2 Antibody, HRP conjugated

Product Code	CSB-PA002333LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P16615
Immunogen	Recombinant Human Sarcoplasmic/endoplasmic reticulum calcium ATPase 2 protein (314-756AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Sarcoplasmic/endoplasmic reticulum calcium ATPase 2 (SERCA2) (SR Ca(2+)-ATPase 2) (EC 3.6.3.8) (Calcium pump 2) (Calcium-transporting ATPase sarcoplasmic reticulum type, slow twitch skeletal muscle isoform) (Endoplasmic reticulum class 1/2 Ca(2+) ATPase), ATP2A2, ATP2B
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	ATP2A2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.