



PAX3 Antibody, FITC conjugated

Product Code	CSB-PA017489HC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P23760
Immunogen	Recombinant Human Paired box protein Pax-3 protein (90-184AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Paired box protein Pax-3 (HuP2), PAX3, HUP2
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	PAX3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.