





JUP Antibody

Product Code	CSB-PA12839A0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P14923
Immunogen	Recombinant Human Junction plakoglobin protein (1-131AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:500-1:5000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Junction plakoglobin (Catenin gamma) (Desmoplakin III) (Desmoplakin-3), JUP, CTNNG DP3
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	JUP

Image



Western blot
 All lanes: JUP antibody at 2µg/ml + A431 whole cell lysate
 Secondary
 Goat polyclonal to rabbit IgG at 1/10000 dilution
 Predicted band size: 82 kDa
 Observed band size: 82 kDa



Immunohistochemistry of paraffin-embedded human breast cancer using CSB-PA12839A0Rb at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.