



Sumo tag Monoclonal Antibody

Product Code	CSB-MA000132M0m
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Immunogen	sumo tag protein
Raised In	mouse
Species Reactivity	N/A
Tested Applications	ELISA,WB;Recommended dilution:WB:1:500-1:5000
Form	liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%,Protein G purified
Isotype	IgG1
Clonality	monoclonal
Product Type	Tag Control Antibody
Clone No.	1A52F8

Image



WB:SUMO-tagged fusion protein(20ng) was subjected to SDS-PAGE followed by Western Blot with CSB-MA000132M0m at dilution of Lane 1:1000, Lane 2:2000, Lane 3:4000 Lane 4:8000, Lane 5:16000, Lane 6:32000 Lane 7:64000, Lane 8:128000, Lane 9:256000
 Secondary
 Goat polyclonal to Mouse IgG at 1/5000 dilution
 Predicted band size:55kd
 Observed band size:55kd

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.