



GCC1 Antibody, Biotin conjugated

Product Code	CSB-PA846598LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96CN9
Immunogen	Recombinant Human GRIP and coiled-coil domain-containing protein 1 protein (4-775AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	GRIP and coiled-coil domain-containing protein 1 (Golgi coiled-coil protein 1), GCC1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	GCC1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.