



PDZD3 Antibody, HRP conjugated

Product Code	CSB-PA774807LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q86UT5
Immunogen	Recombinant Human Na(+)/H(+) exchange regulatory cofactor NHE-RF4 protein (42-291AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Na(+)/H(+) exchange regulatory cofactor NHE-RF4 (NHERF-4) (Intestinal and kidney-enriched PDZ protein) (Natrium-phosphate cotransporter IIa C-terminal-associated protein 2) (Na/Pi cotransporter C-terminal-associated protein 2) (NaPi-Cap2) (PDZ domain-containing protein 2) (PDZ domain-containing protein 3) (Sodium-hydrogen exchanger regulatory factor 4), PDZD3, IKEPP NHERF4 PDZK2
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	PDZD3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.