




SPANXD Antibody

Product Code	CSB-PA819444LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BXN6
Immunogen	Recombinant Human Sperm protein associated with the nucleus on the X chromosome D protein (1-97AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Sperm protein associated with the nucleus on the X chromosome D (Cancer/testis antigen 11.4) (CT11.4) (Nuclear-associated protein SPAN-Xd) (SPANX-D) (SPANX family member D), SPANXD, SPANXE
Immunogen Species	Homo sapiens (Human)
Research Area	Cancer
Target Names	SPANXD

Image


Immunohistochemistry of paraffin-embedded human testis tissue using CSB-PA819444LA01HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.