



SPANXD Antibody, HRP conjugated

| | |
|----------------------------|--|
| Product Code | CSB-PA819444LB01HU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | Q9BXN6 |
| Immunogen | Recombinant Human Sperm protein associated with the nucleus on the X chromosome D protein (1-97AA) |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Tested Applications | ELISA |
| Form | Liquid |
| Conjugate | HRP |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Sperm protein associated with the nucleus on the X chromosome D (Cancer/testis antigen 11.4) (CT11.4) (Nuclear-associated protein SPAN-Xd) (SPANX-D) (SPANX family member D), SPANXD, SPANXE |
| Immunogen Species | Homo sapiens (Human) |
| Research Area | Cancer |
| Target Names | SPANXD |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |