



CASP5 Antibody

Product Code	CSB-PA004550LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P51878
Immunogen	Recombinant Human Caspase-5 protein (137-327AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Caspase-5 (CASP-5) (EC 3.4.22.58) (ICE(rel)-III) (Protease ICH-3) (Protease TY) [Cleaved into: Caspase-5 subunit p20; Caspase-5 subunit p10], CASP5, ICH3
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	CASP5
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.