



NME6 Antibody, FITC conjugated

Product Code	CSB-PA015891LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O75414
Immunogen	Recombinant Human Nucleoside diphosphate kinase 6 protein (1-186AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Nucleoside diphosphate kinase 6 (NDK 6) (NDP kinase 6) (EC 2.7.4.6) (Inhibitor of p53-induced apoptosis-alpha) (IPIA-alpha) (nm23-H6), NME6
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	NME6
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.