




GFP Antibody

Product Code	CSB-PA557657
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P42212
Immunogen	Fusion protein of Aequorea Victoria GFP
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA,WB;ELISA:1:2000-1:10000,WB:1:1000-1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Purification Method	Antigen affinity purification
Isotype	IgG
Alias	Green fluorescent protein
Immunogen Species	Aequorea victoria (Jellyfish)
Target Names	GFP

Image



Gel: 8%SDS-PAGE, Lysate: 1,5,10,20 µg, Lane: GFP transfected 293T cells, Primary antibody: CSB-PA557657(GFP Antibody) at dilution 1/800, Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution, Exposure time: 30 seconds

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.