




STAT6 Antibody

Product Code	CSB-PA138665
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P42226
Immunogen	Fusion protein of Human STAT6
Raised In	Rabbit
Species Reactivity	Human,Mouse
Tested Applications	ELISA,WB;ELISA:1:2000-1:5000,WB:1:500-1:2000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Purification Method	Antigen affinity purification
Isotype	IgG
Alias	signal transducer and activator of transcription 6, interleukin-4 induced
Immunogen Species	Homo sapiens (Human)
Target Names	STAT6

Image



Gel: 8%SDS-PAGE, Lysate: 40 µg, Lane: Raji cells, Primary antibody: CSB-PA138665(STAT6 Antibody) at dilution 1/450, Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution, Exposure time: 5 seconds

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.