




GST Tag Antibody

Product Code	CSB-PA255432
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Immunogen	Full length fusion protein
Raised In	Rabbit
Species Reactivity	GST Tag
Tested Applications	ELISA,WB;ELISA:1:1000-1:10000,WB:1:1000-1:8000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Purification Method	Antigen affinity purification
Isotype	IgG
Alias	Glutathione S-Transferase

Image



Gel: 10%SDS-PAGE, Lysate: 2 µg, Lane: Recombinant GST tag protein, Primary antibody: CSB-PA255432(GST Tag Antibody) at dilution 1/1000, Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution, Exposure time: 15 seconds

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.