




S100A4 Antibody

Product Code	CSB-PA286071
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P26447
Immunogen	Fusion protein of Human S100A4
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,IHC;ELISA:1:500-1:5000,IHC:1:5-1:20
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Purification Method	Antigen affinity purification
Isotype	IgG
Alias	Protein S100-A4
Immunogen Species	Homo sapiens (Human)
Target Names	S100A4

Image



The image on the left is immunohistochemistry of paraffin-embedded Human stomach cancer tissue using CSB-PA286071(S100A4 Antibody) at dilution 1/100, on the right is treated with fusion protein. (Original magnification: ×200)