





BOLL Antibody

Product Code	CSB-PA177370
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8N9W6
Immunogen	Fusion protein of Human BOLL
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, IHC; ELISA: 1:500-1:5000, IHC: 1:5-1:20
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Purification Method	Antigen affinity purification
Isotype	IgG
Alias	Protein boule-like
Immunogen Species	Homo sapiens (Human)
Target Names	BOLL

Image



The image is immunohistochemistry of paraffin-embedded Human lung cancer tissue using CSB-PA177370(BOLL Antibody) at dilution 1/10. (Original magnification: x200)



The image is immunohistochemistry of paraffin-embedded Human ovarian cancer tissue using CSB-PA177370(BOLL Antibody) at dilution 1/10. (Original magnification: x200)

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.