




RBMY1A1 Antibody

Product Code	CSB-PA289764
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P0DJJ3
Immunogen	Fusion protein of Human RBMY1A1
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA,WB;ELISA:1:2000-1:5000,WB:1:500-1:2000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Purification Method	Antigen affinity purification
Isotype	IgG
Alias	RNA binding motif protein, Y-linked, family 1, member A1
Immunogen Species	Homo sapiens (Human)
Target Names	RBMY1A1

Image



Gel: 6%SDS-PAGE, Lysate: 40 µg, Lane 1-4: A549 cells, Hela cells, PC3 cells, HepG2 cells, Primary antibody: CSB-PA289764(RBMY1A1 Antibody) at dilution 1/360, Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution, Exposure time: 30 seconds

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.