




RBMY1A1 Antibody

| | |
|----------------------------|---|
| Product Code | CSB-PA586044 |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | P0DJD3 |
| Immunogen | Fusion protein of Human RBMY1A1 |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Tested Applications | ELISA, WB; ELISA: 1:2000-1:5000, WB: 1:500-1:2000 |
| Form | Liquid |
| Conjugate | Non-conjugated |
| Storage Buffer | -20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol |
| Purification Method | Antigen affinity purification |
| Isotype | IgG |
| Alias | RNA binding motif protein, Y-linked, family 1, member A1 |
| Immunogen Species | Homo sapiens (Human) |
| Target Names | RBMY1A1 |

Image



Gel: 6%SDS-PAGE, Lysate: 40 µg, Lane 1-4: A549 cells, Hela cells, PC3 cells, HepG2 cells, Primary antibody: CSB-PA586044(RBMY1A1 Antibody) at dilution 1/350, Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution, Exposure time: 30 seconds