




ACTG1 Antibody

Product Code	CSB-PA218620
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P63261
Immunogen	Fusion protein of Human ACTG1
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB;ELISA:1:2000-1:5000,WB:1:500-1:2000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Purification Method	Antigen affinity purification
Isotype	IgG
Alias	actin, gamma 1
Immunogen Species	Homo sapiens (Human)
Target Names	ACTG1

Image



Gel: 10+12%SDS-PAGE, Lysate: 40 µg, Lane: Human lymphoma tissue, Primary antibody: CSB-PA218620(ACTG1 Antibody) at dilution 1/350, Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution, Exposure time: 1 second

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.