





ADD1 Antibody

Product Code	CSB-PA985327
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P35611
Immunogen	Fusion protein of Human ADD1
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,IHC;ELISA:1:500-1:5000,IHC:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Purification Method	Antigen affinity purification
Isotype	IgG
Alias	Alpha-adducin
Immunogen Species	Homo sapiens (Human)
Target Names	ADD1

Image



The image is immunohistochemistry of paraffin-embedded Human cervical cancer tissue using CSB-PA985327(ADD1 Antibody) at dilution 1/50. (Original magnification: x200)



The image is immunohistochemistry of paraffin-embedded Human lung cancer tissue using CSB-PA985327(ADD1 Antibody) at dilution 1/50. (Original magnification: x200)

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.