





AFP Antibody

| | |
|----------------------------|--|
| Product Code | CSB-PA230953 |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | P02771 |
| Immunogen | Fusion protein of Human AFP |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Tested Applications | ELISA,WB,IHC;ELISA:1:500-1:10000,WB:1:1000-1:5000,IHC:1:50-1:200 |
| Form | Liquid |
| Conjugate | Non-conjugated |
| Storage Buffer | 0.1mol/L NaCl, 15mmol/L NaN ₃ , PH7.2 |
| Purification Method | Antigen affinity purification |
| Isotype | IgG |
| Alias | Alpha-fetoprotein |
| Immunogen Species | Homo sapiens (Human) |
| Target Names | AFP |


Image



The image is immunohistochemistry of paraffin-embedded Human thyroid cancer tissue using CSB-PA230953(AFP Antibody) at dilution 1/40. (Original magnification: X200)



The image is immunohistochemistry of paraffin-embedded Human liver cancer tissue using CSB-PA230953(AFP Antibody) at dilution 1/40. (Original magnification: X200)



Gel: 8%SDS-PAGE, Lysate: 40 μ g, Lane: HepG2 cells, Primary antibody: CSB-PA230953(AFP Antibody) at dilution 1/1000, Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution, Exposure time: 1 minute

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.