




AGR2 Antibody

Product Code	CSB-PA960585
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O95994
Immunogen	Fusion protein of Human AGR2
Raised In	Rabbit
Species Reactivity	Human,Mouse
Tested Applications	ELISA,WB,IHC;ELISA:1:1000-1:10000,WB:1:1000-1:5000,IHC:1:10-1:50
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Purification Method	Antigen affinity purification
Isotype	IgG
Alias	Anterior gradient protein 2 homolog
Immunogen Species	Homo sapiens (Human)
Target Names	AGR2

Image



The image on the left is immunohistochemistry of paraffin-embedded Human renal cancer tissue using CSB-PA960585(AGR2 Antibody) at dilution 1/10, on the right is treated with fusion protein. (Original magnification: x200)



Gel: 10%SDS-PAGE, Lysate: 40 μg, Lane: Mouse Large intestine tissue, Primary antibody: CSB-PA960585(AGR2 Antibody) at dilution 1/1000, Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution, Exposure time: 20 seconds

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.