





AKR1B1 Antibody

Product Code	CSB-PA078927
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P15121
Immunogen	Fusion protein of Human AKR1B1
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IHC;ELISA:1:500-1:5000,WB:1:500-1:2000,IHC:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Purification Method	Antigen affinity purification
Isotype	IgG
Alias	aldo-keto reductase family 1, member B1 (aldose reductase)
Immunogen Species	Homo sapiens (Human)
Target Names	AKR1B1

Image



The image on the left is immunohistochemistry of paraffin-embedded Human lung cancer tissue using CSB-PA078927(AKR1B1 Antibody) at dilution 1/40, on the right is treated with fusion protein. (Original magnification: x200)



Gel: 10%SDS-PAGE, Lysate: 40 µg, Lane: Hela cells, Primary antibody: CSB-PA078927(AKR1B1 Antibody) at dilution 1/340, Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution, Exposure time: 1 minute