





ALOX5 Antibody

Product Code	CSB-PA043268
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P09917
Immunogen	Fusion protein of Human ALOX5
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,IHC;ELISA:1:500-1:5000,IHC:1:50-1:150
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Purification Method	Antigen affinity purification
Isotype	IgG
Alias	Arachidonate 5-lipoxygenase
Immunogen Species	Homo sapiens (Human)
Target Names	ALOX5

Image



The image is immunohistochemistry of paraffin-embedded Human colon cancer tissue using CSB-PA043268(ALOX5 Antibody) at dilution 1/50. (Original magnification: ×200)



The image is immunohistochemistry of paraffin-embedded Human thyroid cancer tissue using CSB-PA043268(ALOX5 Antibody) at dilution 1/50. (Original magnification: ×200)