





ANXA7 Antibody

Product Code	CSB-PA118689
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P20073
Immunogen	Fusion protein of Human ANXA7
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, IHC; ELISA: 1:1000-1:8000, IHC: 1:150-1:300
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Purification Method	Antigen affinity purification
Isotype	IgG
Alias	Annexin A7
Immunogen Species	Homo sapiens (Human)
Target Names	ANXA7

Image



The image is immunohistochemistry of paraffin-embedded Human cervical cancer tissue using CSB-PA118689(ANXA7 Antibody) at dilution 1/150. (Original magnification: ×200)



The image is immunohistochemistry of paraffin-embedded Human liver cancer tissue using CSB-PA118689(ANXA7 Antibody) at dilution 1/150. (Original magnification: ×200)

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.