





# CASP3 Antibody

<b>Product Code</b>	CSB-PA984107
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	P42574
<b>Immunogen</b>	Fusion protein of Human CASP3
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human,Mouse,Rat
<b>Tested Applications</b>	ELISA,WB,IHC;ELISA:1:1000-1:2000,WB:1:200-1:1000,IHC:1:50-1:200
<b>Form</b>	Liquid
<b>Conjugate</b>	Non-conjugated
<b>Storage Buffer</b>	-20°C, pH7.4 PBS, 0.05% NaN <sub>3</sub> , 40% Glycerol
<b>Purification Method</b>	Antigen affinity purification
<b>Isotype</b>	IgG
<b>Alias</b>	caspase 3, apoptosis-related cysteine peptidase
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Target Names</b>	CASP3


## Image



The image is immunohistochemistry of paraffin-embedded Human gastric cancer tissue using CSB-PA984107(CASP3 Antibody) at dilution 1/30. (Original magnification: x200)



The image is immunohistochemistry of paraffin-embedded Human thyroid cancer tissue using CSB-PA984107(CASP3 Antibody) at dilution 1/30. (Original magnification: x200)



Gel: 10%SDS-PAGE,Lysate: 40  $\mu$ g,Lane 1-2:  
Jurkat cells, 293T cells,Primary antibody: CSB-  
PA984107(CASP3 Antibody) at dilution 1/300  
dilution,Secondary antibody: Goat anti rabbit IgG  
at 1/8000 dilution,Exposure time: 10 seconds

**Usage**

For Research Use Only. Not for use in diagnostic or therapeutic procedures.