





CASP2 Antibody

Product Code	CSB-PA076198
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P42575
Immunogen	Fusion protein of Human CASP2
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IHC;ELISA:1:1000-1:10000,WB:1:1000-1:2000,IHC:1:100-1:300
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Purification Method	Antigen affinity purification
Isotype	IgG
Alias	Caspase-2
Immunogen Species	Homo sapiens (Human)
Target Names	CASP2

Image



The image on the left is immunohistochemistry of paraffin-embedded Human urinary bladder tissue using CSB-PA076198(CASP2 Antibody) at dilution 1/100, on the right is treated with fusion protein. (Original magnification: x200)



Gel: 10%SDS-PAGE, Lysate: 30 µg, Lane 1-2: Human fetal kidney tissue, Human lung tissue, Primary antibody: CSB-PA076198(CASP2 Antibody) at dilution 1/950, Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution, Exposure time: 2 minutes

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.