





CASP7 Antibody

Product Code	CSB-PA505238
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P55210
Immunogen	Fusion protein of Human CASP7
Raised In	Rabbit
Species Reactivity	Human,Mouse
Tested Applications	ELISA,WB,IHC;ELISA:1:1000-1:2000,WB:1:200-1:1000,IHC:1:25-1:100
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Purification Method	Antigen affinity purification
Isotype	IgG
Alias	caspase 7, apoptosis-related cysteine peptidase
Immunogen Species	Homo sapiens (Human)
Target Names	CASP7


Image



The image on the left is immunohistochemistry of paraffin-embedded Human liver cancer tissue using CSB-PA505238(CASP7 Antibody) at dilution 1/20, on the right is treated with fusion protein. (Original magnification: X200)



The image on the left is immunohistochemistry of paraffin-embedded Human lung cancer tissue using CSB-PA505238(CASP7 Antibody) at dilution 1/20, on the right is treated with fusion protein. (Original magnification: X200)



Gel: 8+12%SDS-PAGE, Lysate: 40 μ g, Lane:
 Mouse brain tissue, Primary antibody: CSB-
 PA505238(CASP7 Antibody) at dilution 1/186,
 Secondary antibody: Goat anti rabbit IgG at
 1/8000 dilution, Exposure time: 1 minute