




CASP14 Antibody

Product Code	CSB-PA557650
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P31944
Immunogen	Fusion protein of Human CASP14
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA,WB;ELISA:1:1000-1:2000,WB:1:200-1:1000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Purification Method	Antigen affinity purification
Isotype	IgG
Alias	caspase 14, apoptosis-related cysteine peptidase
Immunogen Species	Homo sapiens (Human)
Target Names	CASP14

Image



Gel: 12%SDS-PAGE, Lysate: 40 µg, Lane 1-3: Jurkat cells, A431 cells, A375 cells, Primary antibody: CSB-PA557650(CASP14 Antibody) at dilution 1/300, Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution, Exposure time: 20 seconds

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.