




TP53 Antibody

Product Code	CSB-PA150453
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P04637
Immunogen	Fusion protein of Human TP53
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB; ELISA: 1:1000-1:5000, WB: 1:500-1:2000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Purification Method	Antigen affinity purification
Isotype	IgG
Alias	tumor protein p53
Immunogen Species	Homo sapiens (Human)
Target Names	TP53

Image



Gel: 12%SDS-PAGE, Lysate: 40 µg, Lane 1-2: Human breast tissue, Human colon cancer tissue, Primary antibody: CSB-PA150453(TP53 Antibody) at dilution 1/300, Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution, Exposure time: 60 seconds