





MMP3 Antibody

Product Code	CSB-PA014205
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P08254
Immunogen	Fusion protein of Human MMP3
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA,IHC;ELISA:1:500-1:5000,IHC:1:25-1:100
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Purification Method	Antigen affinity purification
Isotype	IgG
Alias	Matrix metalloproteinase-3
Immunogen Species	Homo sapiens (Human)
Target Names	MMP3

Image



The image is immunohistochemistry of paraffin-embedded Human breast cancer tissue using CSB-PA014205(MMP3 Antibody) at dilution 1/10. (Original magnification: x200)



The image is immunohistochemistry of paraffin-embedded Human cervical cancer tissue using CSB-PA014205(MMP3 Antibody) at dilution 1/10. (Original magnification: x200)

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.