




TLR7 Antibody

Product Code	CSB-PA792911
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NYK1
Immunogen	Fusion protein of Human TLR7
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, IHC; ELISA: 1:500-1:5000, IHC: 1:10-1:50
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Purification Method	Antigen affinity purification
Isotype	IgG
Alias	Toll-like receptor 7
Immunogen Species	Homo sapiens (Human)
Target Names	TLR7

Image



The image is immunohistochemistry of paraffin-embedded Human tonsil tissue using CSB-PA792911 (TLR7 Antibody) at dilution 1/10. (Original magnification: 200)

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.