




COX5B Antibody

Product Code	CSB-PA573665
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P10606
Immunogen	Fusion protein of Human COX5B
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB;ELISA:1:1000-1:5000,WB:1:200-1:1000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Purification Method	Antigen affinity purification
Isotype	IgG
Alias	cytochrome c oxidase subunit Vb
Immunogen Species	Homo sapiens (Human)
Target Names	COX5B

Image



Gel: 10+12%SDS-PAGE, Lysate: 40 µg, Lane 1-2: Mouse heart tissue, human liver cancer tissue, Primary antibody: CSB-PA573665(COX5B Antibody) at dilution 1/350, Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution, Exposure time: 40 seconds

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.