





COX7B Antibody

Product Code	CSB-PA017662
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P24311
Immunogen	Fusion protein of Human COX7B
Raised In	Rabbit
Species Reactivity	Human,Rat
Tested Applications	ELISA,WB,IHC;ELISA:1:2000-1:5000,WB:1:500-1:2000,IHC:1:100-1:300
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Purification Method	Antigen affinity purification
Isotype	IgG
Alias	Cytochrome c oxidase subunit VIIb
Immunogen Species	Homo sapiens (Human)
Target Names	COX7B


Image



The image on the left is immunohistochemistry of paraffin-embedded Human colon cancer tissue using CSB-PA017662(COX7B Antibody) at dilution 1/70, on the right is treated with fusion protein. (Original magnification: x200)



The image on the left is immunohistochemistry of paraffin-embedded Human gastric cancer tissue using CSB-PA017662(COX7B Antibody) at dilution 1/70, on the right is treated with fusion protein. (Original magnification: x200)



Gel: 10+15%SDS-PAGE, Lysate: 50 μ g, Lane:
293T cells, Primary antibody: CSB-
PA017662(COX7B Antibody) at dilution 1/800,
Secondary antibody: Goat anti rabbit IgG at
1/8000 dilution, Exposure time: 40 seconds

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.