





# COX6B1 Antibody

<b>Product Code</b>	CSB-PA877027
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	P14854
<b>Immunogen</b>	Full length fusion protein
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human,Mouse
<b>Tested Applications</b>	ELISA,WB,IHC;ELISA:1:2000-1:10000,WB:1:1000-1:5000,IHC:1:50-1:200
<b>Form</b>	Liquid
<b>Conjugate</b>	Non-conjugated
<b>Storage Buffer</b>	-20°C, pH7.4 PBS, 0.05% NaN <sub>3</sub> , 40% Glycerol
<b>Purification Method</b>	Antigen affinity purification
<b>Isotype</b>	IgG
<b>Alias</b>	Cytochrome c oxidase subunit VIb polypeptide 1 (ubiquitous)
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Target Names</b>	COX6B1


## Image



The image is immunohistochemistry of paraffin-embedded Human thyroid cancer tissue using CSB-PA877027(COX6B1 Antibody) at dilution 1/50. (Original magnification: x200)



The image is immunohistochemistry of paraffin-embedded Human brain tissue using CSB-PA877027(COX6B1 Antibody) at dilution 1/50. (Original magnification: x200)



Gel: 10%SDS-PAGE,Lysate: 40  $\mu$ g,,Primary antibody: CSB-PA877027(COX6B1 Antibody) at dilution 1/1200 dilution,Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution,Exposure time: 1 minute

## Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.