





APP Antibody

Product Code	CSB-PA891395
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P05067?
Immunogen	Fusion protein of Human APP
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA,IHC;ELISA:1:1000-1:2000,IHC:1:25-1:100
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Purification Method	Antigen affinity purification
Isotype	IgG
Alias	amyloid beta (A4) precursor protein
Immunogen Species	Homo sapiens (Human)
Target Names	APP

Image



The image on the left is immunohistochemistry of paraffin-embedded Human breast cancer tissue using CSB-PA891395(APP Antibody) at dilution 1/25, on the right is treated with fusion protein. (Original magnification: x200)



The image on the left is immunohistochemistry of paraffin-embedded Human thyroid cancer tissue using CSB-PA891395(APP Antibody) at dilution 1/25, on the right is treated with fusion protein. (Original magnification: x200)