





SNCG Antibody

Product Code	CSB-PA580970
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O76070
Immunogen	Fusion protein of Human SNCG
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA,WB,IHC;ELISA:1:2000-1:10000,WB:1:500-1:2000,IHC:1:100-1:300
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Purification Method	Antigen affinity purification
Isotype	IgG
Alias	synuclein, gamma (breast cancer-specific protein 1)
Immunogen Species	Homo sapiens (Human)
Target Names	SNCG


Image



The image is immunohistochemistry of paraffin-embedded Human brain tissue using CSB-PA580970(SNCG Antibody) at dilution 1/60. (Original magnification: x200)



The image is immunohistochemistry of paraffin-embedded Human ovarian cancer tissue using CSB-PA580970(SNCG Antibody) at dilution 1/60. (Original magnification: x200)



Gel: 12%SDS-PAGE, Lysate: 40 μ g, Lane: HT29 cells, Primary antibody: CSB-PA580970(SNCG Antibody) at dilution 1/800, Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution, Exposure time: 20 seconds