





SNCG Antibody

Product Code	CSB-PA830872
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O76070
Immunogen	Fusion protein of Human SNCG
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA,IHC;ELISA:1:2000-1:5000,IHC:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Purification Method	Antigen affinity purification
Isotype	IgG
Alias	synuclein, gamma (breast cancer-specific protein 1)
Immunogen Species	Homo sapiens (Human)
Target Names	SNCG

Image



The image is immunohistochemistry of paraffin-embedded Human thyroid cancer tissue using CSB-PA830872(SNCG Antibody) at dilution 1/50. (Original magnification: x200)



The image is immunohistochemistry of paraffin-embedded Human liver cancer tissue using CSB-PA830872(SNCG Antibody) at dilution 1/50. (Original magnification: x200)

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.