




NFKB1 Antibody

Product Code	CSB-PA259078
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P19838
Immunogen	Fusion protein of Human NFKB1
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,IHC;ELISA:1:500-1:5000,IHC:1:5-1:20
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Purification Method	Antigen affinity purification
Isotype	IgG
Alias	nuclear factor of kappa light polypeptide gene enhancer in B-cells 1
Immunogen Species	Homo sapiens (Human)
Target Names	NFKB1

Image



The image is immunohistochemistry of paraffin-embedded Human thyroid cancer tissue using CSB-PA259078(NFKB1 Antibody) at dilution 1/5. (Original magnification: 200)

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.